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Collagen Gel

Product Introduction

Collagen Gel is a filler obtained from allograft collagen tissue sources that is prepared in the form of lyophilized, cross-linked and irradiated. Due to the main source of this filler, which is mainly obtained from collagen from the dermal layer of human skin, also due to the presence of dense collagen fibers type 1 and 3 in it and the strengthening of these fibers following a specific cross-linking process, it is more durable than other existing collagen gels.

Biocompatibility

This biocompatible product provides adaptable pH with human physiological pH which is suitable for stem cell and fibroblast seeding.

Processing Method & Sterilization

After screening healthy donors with age limitation (more than 13 years old) by evaluating medical history and physical examination, we have quality control tests with microbiologic & serologic tests. Each sample is tested for aerobic and anaerobic bacteria, fungi & HIV antibodies type 1, 2 (anti- HIV1&2), antibody to human T-lymphotropic virus type 1,2 (HTLV Ab1,2), Hepatitis B core antibody (HBcAb), Hepatitis C virus antibody (HCVAb), Hepatitis B surface Antigen (HBsAg), & RPR for syphilis. The confirmation of HIV, HCV and HBV PCR tests is also mandatory. Before cleaning & processing, final validation has been carried out. All procedure is done under the clean conditions of the class B cleanroom at I.T.P Co. After complete processing and packaging procedure, we also use Gamma irradiation for final sterilization.

Indications for Use

As a filler to fill medium and deep lines of the face (cosmetic surgery), correct soft-tissue defects, depressed scars and acne, replacement of integumental tissue lost through atrophy, restore volume of facial components (maxillofacial surgery, ear, nose and throat surgery), use in vocal cord paralysis, incontinency, to fill the laugh line (Nasal Area).

Contraindications

Collagen gel should not be used in existence of infection/skin disease at the surgical site pre-injection of page gel in place, in the forehead, glabellar areas or fine wrinkles around the eyes.

Side Effects

As with any surgical procedure, adverse reactions are possible and may include, but are not limited swelling, bruising, pain or irritation during the immediate post-operative period. Skin redness, discoloration and herpetic lesions have been reported on rare occasions. These symptoms typically resolve without sequelae.

Presentation and Packaging

Collagen Gel in syringe in the double film pocket that placed in film and medical grade paper.

Warning and Precautions

Do not use this product if:

- ✓ The expiration date on the package is passed

- ✓ Any of the package or product appears to be missed, discolored or damaged.
- ✓ The product label is damaged and missing.

Caution

Please use the contents of any package for only one patient. Please avoid resterilization the product.

Please report any adverse effect attributable to this product to us as soon as possible.

Contact Number: +98 9024264707

Storage

Store this product at refrigerator (2-8 °C). Products have 6 months shelf life.

Instructions for Use

- ✓ Peel open outer package (nonsterile).
- ✓ Using sterile technique, transfer contents to a sterile field.
- ✓ Remove protective cap from syringe.
- ✓ Inject product in the site.
- ✓ Discard any unused portion.

Health Insurance Products

- ✓ Implementation of instructions for product preparation and processing in compliance with the requirements of the Food and Drug Administration of Iran and the National Medical Device Directorate.
- ✓ Observance of GMP production principles
- ✓ Obtaining the international standard for medical devices- Quality management systems ISO 13485
- ✓ Receiving the product license and originality from National Medical Device Directorate labels from the Food and Drug Administration of Iran.
- ✓ Supervising the quality control system In-House and supervising the competent institutions of the Ministry of Health and medical education on preparation, production and distribution.

Manufacturer:

Address: Iranian Tissue Product Co.

Danesh 2 St., Phase 2, Pardis Technology Park, 20th km Damavand Road, Tehran, Iran.

Telefax: +9821 68105

E-mail address: info@regen.ir

References

The 14th Edition of American Association of Tissue Banks Standards for tissue banking (AATB)

-European Good Tissue & Cell Practices (Euro GTP)