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Product Introduction

Bone Allograft Chips, Cube, Matchstick, Block and Wedge Mineralized Freeze Dried is processed from cancellous or cortical - cancellous bone. Mineralized bone matrices provide a biocompatible osteoconductive matrix that support new bone formation.

Biocompatibility

This biocompatible product provides adaptable pH with human physiological pH which is suitable for stem cell and osteoprogenitor seeding.

Processing Method & Sterilization

After screening healthy donors with age limitation (6 to 70 years old) by evaluating medical history and physical examination, we have quality control tests with microbiologic & serologic tests. Each sample is tested for aerobic and anaerobic bacteria, fungi & HIV antibodies type 1, 2 (anti- HIV1&2), antibody to human T-lymphotropic virus type 1,2 (HTLV Ab1,2), Hepatitis B core antibody (HBcAb), Hepatitis C virus antibody (HCVAb), Hepatitis B surface Antigen (HBsAg), & RPR for syphilis. The confirmation of HIV, HCV and HBV PCR tests is also mandatory. Before cleaning & processing, final validation has been carried out. All procedure is done under the clean conditions of the class C cleanroom at I.T.P Co. After complete processing and packaging procedure, we also use Gamma irradiation for final sterilization.

Indications for Use

It is used for reconstruction or filling bone defects in dentistry, periodontology, oral and maxillofacial surgery, reconstructive surgery.

Contraindications

Do not use product when the necrotic tissue, active or latent infection is in or around the surgical site. Do not use this product on unknown defects too.

Side Effects

As with any surgical procedure, adverse reactions are possible and may include, but are not limited to infection, inflammation, and rejection.

Presentation and Packaging

MBA Bone Substitute packed in sterilized single use package and has been presented for single patient use. The product has been packed in double pouches, the first pouches from double film that placed in film and medical grade paper. The product has been sterilized through gamma irradiation method.

Warning and Precautions

Do not use this product if:

- ✓ The expiration date on the package is passed
- ✓ Any of the package or product appears to be missed, discolored or damaged.
- ✓ The product label is damaged and missing.

Caution

Please use the contents of any package for only one patient. Please avoid resterilization the product.

Please report any adverse effect attributable to this product to us as soon as possible.

Contact Number: +98 9024264707

Storage

Store this product at ambient temperature (15-30 °C). No refrigeration or freezing is needed. products have 5 years shelf life.

Preoperative preparation

Preparation of the bone graft bed is important for graft in corporation and bone formation. Sterile technique must be maintained to minimize the risk of postoperative complications. The amount of graft material needed is based on the defect size and the procedure & is determined by the judgment of the clinician.

Instructions for Use

- ✓ Peel open outer package (nonsterile).
- ✓ Using sterile technique, transfer contents to a sterile field.
- ✓ Soak product in a sustainable physiological medium or recipient's blood for about 15-20 minutes (mandatory).
- ✓ Put the combination in the surgical site.
- ✓ Discard any unused portion.

Health Insurance Products

- ✓ Implementation of instructions for product preparation and processing in compliance with the requirements of the Food and Drug Administration of Iran and the National Medical Device Directorate.
- ✓ Observance of GMP production principles.
- ✓ Obtaining the international standard for medical devices- Quality management systems ISO 13485
- ✓ Receiving the product license and originality from National Medical Device Directorate labels from the Food and Drug Administration of Iran.
- ✓ Supervising the quality control system In-House and supervising the competent institutions of the Ministry of Health and medical education on preparation, production and distribution.

Manufacturer:

Address: Iranian Tissue Product Co.

Danesh 2 St., Phase 2, Pardis Technology Park, 20th km Damavand Road, Tehran, Iran.

Telefax: +9821 68105

E-mail address: info@regen.ir

References

- The 14th Edition of American Association of Tissue Banks Standards for tissue banking (AATB)
- European Good Tissue & Cell Practices (Euro GTP)