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Product Introduction

The septum is a cartilage-bone wall that separates the left and right nasal cavities, which removes the required cartilage from the septum and thereby reinforces the tip of the nose. Cartilage allograft is processed from donated human cartilage tissue. We preserve allogeneous septum in ethanol for easier use.

Biocompatibility

This biocompatible product provides adaptable pH with human physiological pH which is suitable for stem cell and Chondrocytes seeding.

Processing Method & Sterilization

After screening healthy donors with age limitation (15 to 60 years old) by evaluating medical history and physical examination, we have quality control tests with microbiologic & serologic tests. Each sample is tested for aerobic and anaerobic bacteria, fungi & HIV antibodies type 1, 2 (anti- HIV1&2), antibody to human T-lymphotropic virus type 1,2 (HTLV Ab1,2), Hepatitis B core antibody (HBcAb), Hepatitis C virus antibody (HCVAb), Hepatitis B surface Antigen (HBsAg), & RPR for syphilis. The confirmation of HIV, HCV and HBV PCR tests is also mandatory. Before cleaning & processing, final validation has been carried out. All procedure is done under clean conditions of class B cleanroom at I.T.P Co.

Indications for Use

It is used for septorhinoplasty surgery, oral and maxillofacial surgery, E.N.T & plastic and reconstructive surgery.

Contraindications

Do not use product when the necrotic tissue, active or latent infection is in or around the surgical site. Do not use this product on unknown defects too.

Side Effects

There is a possibility of infection, inflammation, adverse reaction and rejection.

Presentation and Packaging

The package should be inspected before use to ensure that it is not damaged. Septum stored in sterilized piriform and double pouched film pocket.

Warning and Precautions

Do not use this product if:

- ✓ The expiration date on the package is passed
- ✓ Any of the package or product appears to be missed, discolored or damaged.
- ✓ The product label is damaged and missing.

Caution

Please use the contents of any package for only one patient. Please avoid resterilization the product.

Please report any adverse effect attributable to this product to us as soon as possible.

Contact Number: +98 9024264707

Storage

Store Alcohol preserved septum in refrigerator (2-8 °C). It is necessary to check the volume of the alcohol prior to use. products have 2-year shelf life.

Preoperative preparation

Preparation of the graft bed is important for graft in corporation and new tissue formation. Sterile technique must be maintained to minimize the risk of postoperative complications. The amount of graft material needed is based on the defect size and the procedure & is determined by the judgment of the clinician.

Instructions for Use

- ✓ Peel open outer package (nonsterile).
- ✓ Using sterile technique, transfer contents to a sterile field.
- ✓ Put the combination in the surgical site.
- ✓ Discard any unused portion.

Health Insurance Products

- ✓ Implementation of instructions for product preparation and processing in compliance with the requirements of the Food and Drug Administration of Iran and the National Medical Device Directorate.
- ✓ Observance of GMP production principles
- ✓ Obtaining the international standard for medical devices- Quality management systems ISO 13485
- ✓ Receiving the product license and originality from National Medical Device Directorate labels from the Food and Drug Administration of Iran.
- ✓ Supervising the quality control system In-House and supervising the competent institutions of the Ministry of Health and medical education on preparation, production, and distribution.

Manufacturer:

Address: Iranian Tissue Product Co.

Danesh 2 St., Phase 2, Pardis Technology Park, 20th km Damavand Road, Tehran, Iran.

Telefax: +9821 68105

E-mail address: info@regen.ir

References

- The 14th Edition of American Association of Tissue Banks Standards for tissue banking (AATB)
- European Good Tissue & Cell Practices (Euro GTP)