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Combination of Mineralized and Demineralized Powder and Granule Bone

Product Introduction

Powder and Granule Mineralized and Demineralized Freeze-Dried Bone Allograft is processed from cortical-cancellous bone. Mineralized and demineralized particles and granules are a mixture of 65% MBA and 35% DBM. Mineralized bone matrix provides a biocompatible and osteoconductive environment, and demineralized bone matrix supports osteogenesis and new bone formation.

Biocompatibility

This biocompatible product provides adaptable pH with human physiological pH which is suitable for stem cell and preosteoblast seeding.

Processing Method & Sterilization

After screening healthy donors with age limitation (13 to 70 years old) by evaluating medical history and physical examination, we have quality control tests with microbiologic & serologic tests. Each sample is tested for aerobic and anaerobic bacteria, fungi & HIV antibodies type 1, 2 (HIV Ab I, II), antibody to human T-lymphotropic virus type 1, 2 (HTLV Ab I, II), Hepatitis B core antibody (HBcAb), Hepatitis C virus antibody (HCVAb), Hepatitis B surface Antigen (HBsAg), & RPR for syphilis. The confirmation of HIV, HCV and HBV PCR tests is also mandatory. We have final validation for bone before cleaning & processing. All procedure is done under clean conditions of class C cleanroom at I.T.P Co. After complete processing and packaging procedure, we also use Gamma irradiation for final sterilization.

Indications for Use

It is used for reconstruction or filling bone defects in dentistry, periodontology, oral and maxillofacial surgery, reconstructive surgery.

Contraindications

Do not use product when the necrotic tissue, active or latent infection is in or around the surgical site.

Do not use this product on unknown defects too.

Side Effects

There is a possibility of infection, inflammation, adverse reaction and rejection.

Presentation and Packaging

MBA/DBM Bone Substitute powder and granule are packaged in sterile, single-use packages. The package consists of a glass vial and special paper. This package should be considered to make sure that it has not been defaced.

Storage

Store this product at ambient temperature (15-30 °C). No refrigeration or freezing is needed. products have 5 years shelf life.

Warning and Precautions

Do not use this product if:

- The expiration date on the package is passed
- Any of the package or product appears to be missed, discolored or damaged.

- The product label is damaged and missing.

Caution

- Please use the contents of any package for only one patient.
- Please avoid resterilization the product.
- Please report any adverse effect attributable to this product to us as soon as possible.
- Contact Number: +98 9024264707

Preoperative preparation

Preparation of the bone graft bed is important for graft in corporation and bone formation. Maintaining product sterility is essential to minimize postoperative complications. The amount of graft material needed is based on the defect size and the procedure & is determined by the judgment of the clinician.

Instructions for Use

- Peel open outer package (nonsterile).
- Using sterile technique, transfer contents to a sterile field.
- Soak the product in a suitable physiological serum such as sterile distilled water for 15-20 minutes and then in the recipient's blood if needed (mandatory).
- Put the rehydrated product in the surgical site.
- Discard any unused portion.

Health Insurance Products

- Implementation of instructions for product preparation and processing in compliance with the requirements of the Food and Drug Administration of Iran and the National Medical Device Directorate.
- Observance of GMP production principles.
- Obtaining the international standard for medical devices- Quality management systems ISO 13485
- Receiving the product license and originality from National Medical Device Directorate labels from the Food and Drug Administration of Iran.
- Supervising the quality control system In-House and supervising the competent institutions of the Ministry of Health and medical education on preparation, production and distribution.

Manufacturer:

Address: Iranian Tissue Product Co. Danesh 2 St., Phase 2, Pardis Technology Park, 20th km Damavand Road, Tehran, Iran.

Telefax: +9821 68105

E-mail address: info@regen.ir

References

- The 14th Edition of American Association of Tissue Banks Standards for tissue banking (AATB)
- European Good Tissue & Cell Practices (Euro GTP)